

Anomaly Detection in ECG Signals using a Quantum Autoencoder

Wojciech Czerski, Marcel Biežuński

Project Goal and Approach

- The main goal was to investigate whether quantum models can meaningfully separate normal and anomalous ECG patterns.
- One-class learning - trained on normal ECG signals
- Comparison between classical, quantum-simulated and hardware-based inference

Challenges and Motivation

- High dimensionality and noise in ECG signals
- Classical methods rely on handcrafted features or deep models that require large datasets.
- Quantum computing offers a different representation paradigm, where information is processed in high-dimensional Hilbert spaces, potentially enabling richer feature transformations

Comparison of Approaches

- Classical Approach - Classical autoencoder trained end-to-end on real-valued ECG features, serving as a strong baseline optimized for accuracy and scalability.
- Quantum Simulator - Idealized quantum execution without noise, used to evaluate the representational capacity of the quantum autoencoder independently of hardware imperfections.
- Real IBM device - only evaluation on real IBM Quantum hardware to assess the impact of noise, limited shots, and hardware constraints.

Dataset


- MIT-BIH Arrhythmia Database

(Base Setup)


Train set - 1299 Normal

Test set - 1200 Normal, 964 Anomaly


Single Heartbeat (Time Domain)



Qiskit input data (8 qubits)





ECG Signal Processing Pipeline





Classical Baseline: MLP Autoencoder Architecture



- Feed-forward Neural Network (Multi-Layer Perceptron).
- Structure: 8 (Input) \rightarrow 4 (Latent Space) \rightarrow 8 (Output).
- Objective: Identity mapping ($f(x) \approx x$) with dimensionality reduction.
- Mechanism: Information Bottleneck forces the model to learn robust features.
- Optimization: Further compression to 2 neurons yielded higher accuracy.



Results & Analysis



Results & Analysis



Quantum Architecture Overview

The model consists of two main components:


- **Data encoding circuit**

Classical ECG features are embedded into a quantum state using a feature map.

- **Trainable quantum ansatz**

Ansatz transforms this state before measurement producing probabilistic outputs used for anomaly scoring

Feature Visualisation: 2.2182




Experimental Setup

- The model is trained exclusively on normal ECG segments
- Evaluation is performed on separate test sets containing both normal and anomalous samples
- Identical preprocessing steps are applied across all experiments to ensure fair comparison
- $$L(\theta) = E_{x \sim D_{normal}}[P_{\theta}(anomaly | x)]$$


Implementation details

- **Entire pipeline is implemented in Jupyter Notebook**
 - Ansatz - RealAmplitudes
 - Optimizer - SPSA
 - Feature map - ZFeatureMap
- **Libraries used during the project:**
 - Qiskit 2.3.0 - Quantum Circuit
 - Qiskit algorithms 0.4.0 - Optimizers
 - IBM runtime 0.45.0 - Connecting to the real IBM device
 - Scikit-learn 1.8.0 - Preprocessing, Decomposition




Results and Analysis



Learning Curve Quantum Autoencoder)



Quantum Autoencoder (8 qubits)



Results and Analysis



Scalability Analysis: Dataset Size & Generalization

- **Two scalability axis were tested independently:**

- Sample per patient: 200 → 600
- Number of patients: 1→2→3→5

- **For each configuration we measured:**

- Reconstruction loss
- Detection accuracy
- ROC AUC
- Total number of samples

Scalability Analysis

Hidden Layer Neurons	Number of Patients	Sample Limit / Patient	Summed Samples	Loss	Accuracy	ROC AUC
4	1	200	200	0.038846	81.24	0.9946
4	1	600	600	0.033613	81.56	0.9957
4	2	200	400	0.046569	80.36	0.9917
4	2	600	1200	0.048084	80.64	0.9955
4	3	200	600	0.074402	80.22	0.9818
4	3	600	1800	0.060584	80.96	0.9897
4	5	200	800	0.079987	77.63	0.9762
4	5	600	2062	0.06673	77.87	0.9838

Classical

Quantum

Qubits	Number of Patients	Sample limit/patient	Summed samples	Loss	Accuracy	ROC AUC
8	1	200	200	0.261204	0.5175	0.506906
8	1	600	600	0.200252	0.678373	0.703147
8	2	200	400	0.14366	0.7725	0.838231
8	2	600	1200	0.122035	0.859519	0.88273
8	3	200	499	0.18057	0.75375	0.778138
8	3	600	1299	0.037794	0.766636	0.677064
8	5	200	862	0.273116	0.755	0.769144
8	5	600	2062	0.138173	0.744455	0.729215

Optimizing Classical Architecture

- **Observation:**

- 4 neurons: Slightly better for single-patient data, but performance degrades as patient diversity increases (overfitting to noise).
- 2 neurons: Demonstrated superior robustness. Accuracy remained stable (~81%) event with 5 heterogeneous patients.
- Sample Saturation: Increasing training samples from 200 to 600 per patient yielded diminishing returns (<0.5% accuracy gain). The model converges quickly.
- In general more samples translate into better accuracy and a smaller loss function result.

- **Conclusion:** The ECG signal features are highly compressible. A tighter bottleneck forces the model to learn universal features rather than patient-specific noise. High compression (2 latent neurons) in classical autoencoder is essential for a subject-independent model.

Number of Patients	Latent Neurons = 2	Latent Neurons = 4	2vs4 accuracy difference
1	81.05	81.56	-0.51
2	80.45	80.64	-0.19
3	80.91	80.96	-0.05
5	80.91	77.87	3.04

Optimizing Quantum Autoencoder

- **Observation**

- Gains are visible for small patient sets (1–3 patients), while benefits diminish for higher diversity.
- Increasing the number of qubits from 4 to 8 improve performance and stability

- **Sample saturation**

- Similarly to classical case sample increase yields limited accuracy gains
- This suggests that data diversity dominates over data quantity



Number of Patients	Qubits = 4 (Accuracy %)	Qubits = 8 (Accuracy %)	4 → 8 Accuracy Difference
1	67.01	67.84	+0.83
2	83.83	85.95	+2.12
3	72.37	76.66	+4.29
5	74.49	74.45	−0.04

Real IBM QPU

Experimental setup and constraints:

- Execution on real IBM Quantum hardware (QPU)
- 40 ECG samples executed within a single quantum job
- Inference-only run (training performed on simulator)
- Limited sample size due to hardware access and noise constraints

Results and Analysis



Key Findings

- **Feasibility of Quantum Anomaly Detection**

- Successfully implemented a hybrid pipeline using an 8-qubit autoencoder.
- Demonstrated that quantum circuits can learn structured representations of ECG signals.

- **Anomaly Separation Capability**

- The quantum autoencoder achieved measurable separation between normal and anomalous ECG segments.
- Performance exceeded random baseline, as confirmed by ROC analysis.

- **Scalability Trends**

- Stable behavior observed across different dataset sizes.